



NEWBORN AND INFANT BED MINI/MINI2

Newborn and infant bed MINI/MINI2

is intended for neonatal and paediatric wards of all medical institutions, the MINI2 type is suitable for the care of twins. It is well suited for the ROOMING - IN system, as its design allows it to be placed directly at the mother's bedside.





Mobile, manually easily adjustable in height, four side swivel castors, two equipped with brakes. On the bottom of the chassis there is a storage basket.

With a simple positioning mechanism, the loading surface can be adjusted to both "Trendelenburg", raised or so-called drainage positions.



Optional accessories

Mattress cover fabric or vapour permeable Purplast MINT.

The bed can be supplemented with a heated mattress with electronic temperature control, type **HM-13**, heating module **VM-94**, or high-efficiency LED phototherapy **FL-2010**.

Technical parameters

▪ Dimensions	MINI	500 x 800	mm
	MINI2	650 x 800	mm
▪ Height adjustment		250(600-850)	mm
▪ Positioning settings		+/- 20	°
▪ Weight		cca 14 – 20	kg
▪ Colour version		RAL by agreement	





Alfamedic s.r.o

Tř. 5.května 458/35, 37372 Lišov, Czech Republic

T: 387 994 664

alfamedic@alfamedic.cz, www.alfamedic.cz